

KS-1002A

12 in 1 Skin Care Unit



- Vapozone
- Loop (LED)
- Face Lifting
- Ultrasonic For Face
- Ultrasonic For Eye Area
- Ultrasonic For Body
- Ultrasonic Peeling
- Galvanic
- High Frequency
- Frimator
- Spray
- Vacuum



Vapozone, Loop (LED), Face Lifting, Ultrasonic For Face, Ultrasonic For Eye Are, Ultrasonic For Body, Ultrasonic Peeling, Galvanic, High Frequency, Frimator, Spray and Vacuum functions
It is a stylistic, high-tech unit that meets 12 functions in one device besides its **“LCD Touchscreen”** function.



Vapozone

Vapozone is the first step of skin care and with vapozone feature; you can choose “steam” or “ozone”. Moreover, you can adjust energy level and make aromatherapy with essences. It guarantees a safe use with its auto power off feature.

Loop (LED)

LED light, which does not harm eyes while operating, enables an easy use with its 1/3 magnifying feature, internal stretches and stylish design.


Face Lifting

With KS-1002A Face Lifting headpiece, mature, aging and oily skins can be firmed. Electrode (cathode) and alkaline is applied to the skin. With this treatment nerves are activated, blood circulation is increased, lymph structure is smoothened and skin elasticity is improved. Excessive oil in the skin is reduced and skin excretion is balanced. At the same time, it removes heavy metal ions and residuals in the deeper layers of the skin. It also has a detoxifying feature.

Ultrasonic treatment, which has mechanic, heating and chemical effect, can gently massage to tissues, heat the layer of the tissues, stimulate cells and blood circulation, smoothen tissues, create chemical reaction, strengthen metabolism, affect various enzymes, change Ph value and bio activated substances.

Ultrasonic For Face

When used with solutions that treat different type of acnes (including comedo), the headpiece of the ultrasonic for face can treat acnes. You can get a noticeable result in reducing pigmentation problems and various spots when you use the headpiece with decoloring products.



Ultrasonic For Eye Area

With the headpiece of ultrasonic for eye area, it is possible to treat bag under eye and dark circles. The mechanical effect of the ultrasonic thins the fats and resolves them. Moreover, it stimulates blood circulation and disaggregates water and fat deposits.

Ultrasonic For Body

With mechanic massage effect of the headpiece, it enables a stimulation of blood circulation, resolving fats that cause cellulite and firming of the tissues.

Ultrasonic Peeling

The headpiece of KS-1002A Ultrasonic Peeling enables a deep and painless treatment for acnes and especially for skins with problems. It also helps reducing fine lines and wrinkles as well as smoothing the skin.

Galvanic

With galvanic headpiece you can make treatments for dry skins lack moisture, mature skins that lack nutrition and oily skins. Positron (anode) applies acid combination to the skin and it relaxes nerves, reduces blood flow and redness, improves lymph structure and cell permeability, shrinks the pores and enables the penetration of nutritions that are hard to penetrate. It gives the skin an intensive care.

High Frequency

It stimulates the blood circulation as well as giving the skin elasticity and saturation. It strengthens the cell metabolism. Thanks to its heating effect on skin, it improves cell permeability and penetration potential of the skin. As a result of these features, skin care products penetrate more easily.

Frimatör

With its brush set designed for different treatment types, it is possible to operate in big and small areas of the body like hand, foot, face etc. according to the needs of the patient. Rotation side of the frimator, energy level and treatment time can be adjusted on LCD screen.

Spray

It uses the positive air stream that is produced by the vacuum pump. It sprays treatment liquids (i.e. tonic, solution).

Vacuum

It uses negative pressure air stream of the pump and cleans black spots, inflamed pimples, oil and dirt in the pores with its vacuum function.





KS-1002A Technical Features

| | |
|---------------|--|
| Product Name | Skin Care Unit |
| Product Model | KS-1002A |
| Functions | <ul style="list-style-type: none"> • Vapozone • Loop (LED) • Face Lifting • Ultrasonic For Face • Ultrasonic For Eye Area • Ultrasonic For Body • Ultrasonic Peeling • Galvanic • High Frequency • Frimator • Spray • Vacuum |
| Voltage | AC 220-240V ~ |
| Outler Power | 800 W |